




Confocal Raman Imaging and Correlative Techniques in Life Science



Confocal Raman Imaging in Life Science - Living Cells, Plants, Bacteria & Tissues

With confocal Raman imaging the molecules of a sample can be chemically identified and their distribution can be imaged three-dimensionally. These benefits gain more and more recognition in biological and medical research and confocal Raman microscopy (CRM) becomes a frequently used methods for the clarification of crucial questions in life sciences. Thereby diverse samples are investigated on a regular basis: CRM is applied for measurements in liquids and life cell imaging, to solid samples and soft tissues. By various examples the application note explains the possible use of confocal Raman microscopy and its correlative techniques.


Working Principles

Confocal Raman Imaging


The Raman effect is based on light interacting with the chemical bonds of a sample. Due to vibrations in the chemical bonds this interaction causes a specific energy shift in the back scattered light which appear in a unique Raman spectrum through which the molecular components of a sample can be detected. The confocal Raman Imaging technique combines Raman spectroscopy with a confocal microscope. Thus the spatial distribution of the chemical components within the sample are imaged. High-resolution confocal Raman microscopes acquire the information of a complete Raman spectrum at every image pixel and achieve a lateral resolution at the diffraction limit (circa $\lambda/2$ of the excitation wavelength). A confocal microscope setup is furthermore characterized by an excellent depth resolution and features the generation of 3D Raman images and depth profiles.

Confocal Raman microscopy can be


Raman Spectroscopy



Atomic Force Microscopy



Scanning Near-field Optical Microscopy



coupled with correlative microscopy techniques such as fluorescence microscopy, electron microscopy and atomic force microscopy.

Atomic Force Microscopy (AFM)

In Atomic Force Microscopy the sample is scanned under the tip using a piezo-driven scanning-stage and the topography is displayed as an image. Atomic Force Microscopy provides spatial information parallel and perpendicular to the surface with resolution in the nm range. In addition to topographic high-resolution information, local properties such as adhesion and stiffness can be investigated by analyzing the tip-sample interaction forces. In combination with confocal Raman imaging the topographic information can be linked to molecular information.

Scanning Near-field Optical Microscopy (SNOM)

In Scanning Near-field Optical Mi-

croscopy is a high-resolution imaging technique with an optical resolution in the range of 60 – 100 nm. Therefore the excitation laser light is focused through an aperture with a diameter smaller than the excitation wavelength, resulting in an evanescent field (or near-field) on the far side of the aperture. When the sample is scanned at a small distance below the cantilever tip the optical resolution is limited only by the diameter of the aperture. The transmitted or emitted light is then detected point by point and line by line in order to generate an optical image. The aperture itself is located at the apex of a hollow pyramid on the micro-fabricated WITec SNOM cantilever.

Scanning Electron Microscopy (SEM)

SEM is a high-resolution imaging technique that can determine a sample's topography, morphology, and surface structure. The images are

acquired with nano- or micrometer resolutions and magnification ranges of 10 – 10.000x. Due to the large depth of field, it is also possible to generate a 3D appearance of the surface. SEM is a non-destructive method which allows combinations with other imaging techniques such as confocal Raman imaging. This technique is known as RISE (Raman Imaging and Scanning Electron) Microscopy. Through RISE Microscopy, ultra-structural surface properties can be linked to molecular compound information. This unique combination opens up new dimensions for more comprehensive sample characterization.


Raman microscopy of living cells

To study living cells in their physiological surroundings without damaging


them is a highly sought after capability in life science. Confocal Raman imaging is a nondestructive method with the capability to identify chemical components inside a cell. In this experiment, epithelial rat cells were investigated with the Spectral Imaging Mode of a WITec alpha300 Raman microscopy. A spectrum at every pixel was taken (Scan Range: 40 x 40 μm^2 , 100 x 100 pixel, 10,000 spectra) using a 60x Nikon water objective (NA=1.0). The sample was excited with a 10 mW power 532 nm frequency doubled Nd:YAG laser. Using the integrated video camera, a suitable cell was found (A). The Raman image in (B) shows the integral intensity of the CH-stretching band as indicated by the blue rectangle in (C) which shows one of the 10,000 spectra, acquired with an integration time of 100 ms. The proteins

and lipids can then spectroscopically identified by their particular Raman bands.


In order to optimize the signal to noise ratio of the Raman images, a software tool performs a “fit procedure”: From a small region of interest, a “basic spectrum” is generated by averaging all spectra in this particular region. This basic spectrum can be fitted to the measured spectra of the complete measurement, resulting in an optimized image. The three spectra in (D) represent basic spectra of different regions of the cell with their corresponding images. The blue spectrum corresponds to the mitochondria. The green spectrum is the generated basic spectrum of the ER region (endoplasmatic reticulum). The basic spectrum of the nucleoli region is shown in red.





(A) Video image.



(B) Intensity in the CH stretching band (2800 cm^{-1} - 3000 cm^{-1}).




(C) One of the 10,000 spectra recorded.



(D) Determination of the basis spectra with corresponding images. The colors of the spectra correspond to the colors of the image frames.


Samples courtesy of Dr. Angelika Rück, ILM, Ulm, Germany

To distinguish the components of the cell, the three images generated by the fit procedure were illustrated with different colors and combined in one image (E). The blue features represent the mitochondria and the red areas show the nucleoli inside the nucleolus. The endoplasmic reticulum and other different parts of the cell are clearly visible. Even the nuclear membrane is well-defined. Thus, using confocal Raman imaging, a “color-labelled” image can be created without dyeing the sample.



(E) Color-coded image of the different regions of the cell.


AFM measurement in liquids



The image shows a large area topography AFM scan of a cell culture in fluid.
 Scan range: 250 x 100 µm²
 Maximum measured height: 2.5 µm.

Chromosomal dissection with AFM

AFM allows imaging of DNA in ambient as well as physiological conditions, and several experiments have also demonstrated the capability of the AFM to manipulate biological samples. Human metaphase chromosomes were dissected in the AFM mode of the alpha300 microscope. The chromatid was first cut using the WITec DaVinci Nanolithography package and then imaged with AFM. The topography image shows the chromosome with the cut (3 x 3 µm²) and the corresponding zoom-in of the cut region.



Chromosome after dissection and corresponding zoom-in of the cut.

Raman imaging on fluorescent cells


Fluorescence is a well-established and common method for the investigation of biological samples. For fluorescence experiments the expression of certain proteins is coupled to the expression of fluorescent proteins or fluorophore coupled antibodies against certain proteins are applied.

However, Raman microscopy allows the investigation of the molecular and chemical components of a sample. Thus all components of the cells can be investigated and not just those once which are marked by fluorescence. Nevertheless it is sometimes challenging to identify the area of interest in Raman microscopy. In this case, fluorescence helps to enable the


navigation to the specific region prior to the Raman analysis. Since fluorescence emission interferes with the Raman signal, it is important to select for fluorophores that are compatible with the Raman investigations.

Fixed cells on a microscope slide were studied with the alpha300 Fluorescence-Raman microscope. Cells expressing the Green Fluorescence Protein (GFP) in the cytoplasm were identified in the fluorescence mode (figure A) using a Nikon 40x (NA = 0.6) air objective. Then the same cell was imaged with the integrated video camera (figure B) before imaging the cell in confocal Raman imaging mode with a ZEISS 63x (NA = 1.25) oil immersion objective. The sample was excited with a 532 nm diode laser. A


Raman spectrum at every pixel was taken (scan range: $43 \times 43 \mu\text{m}^2$, 120×120 pixels, 14,000 spectra). With the color-coded Raman image the cell body, the cytoplasm, the nucleus and the nucleoli can be clearly distinguished (figure C). In further analysis, the Raman spectra obtained from the Raman measurement shown in figure C were analyzed by cluster analysis. Therefore the Raman spectra are automatically compared and similar spectra are grouped into several clusters. From this analysis a reduced and simplified Raman image can be displayed (figure D) in which proteins are shown in blue, the nucleus is shown in red, the nucleoli are shown in green and lipids are shown in yellow. The corresponding Raman spectra are shown in figure E.





(A) Fluorescence image of the GFP expressing cells, excitation from above



(B) Video image from above



(C) Color-coded Raman image of the fluorescent cell.






(D) Color-coded cluster analysis image

(E) Corresponding Raman spectra of the four clusters.

Image courtesy of Dr. Claudia Salfli-Happ, ILM, Ulm, Germany





Investigating the intracellular components

Confocal Raman microscopy was used to image intracellular components of macrophages. Macrophages take up low density lipoproteins (LDL) for recycling and removal. If their export capabilities are overextended, macrophages store lipids and develop into foam cells. These foam cells contribute to cardiovascular diseases such as arteriosclerosis. The subcellular lipid distribution of an in vitro THP-1 monocytes cell model was analyzed. The monocytes were differentiated into macrophages, incubated with deuterium labeled lipids and examined with the confocal Raman microscope alpha300 R with a 60x/NA=1.0 water immersion objective, an excitation wavelength of 488 nm, and 5mW laser power at the sample. For the lipid uptake quantification of each cell, the CD (2050-2275 cm^{-1}) to CH (2800 – 3020 cm^{-1}) ratio of the Raman scattering intensities were evaluated. Thereby the CD-stretching vibrations were used as marker for intracellular lipids, while the CH-stretching vibrations reflects the general density of the cells composition. A complete Raman spectrum was acquired at every



(A) Raman image of a macrophage cell incubated with 400 μM of oleic acid for 3 hours, generated from the CH-stretching intensities. The image was recorded with 488 nm excitation using a 60x/NA=1.0 water immersion objective at a step size of 0.5 μm . (B) shows a Raman image reconstructed using a spectral unmixing algorithm, which decomposes the data set of the image into the most dissimilar spectral components. The associated spectral information is plotted in (C). Clearly visible are the Raman signals that originate from the CH-stretching vibrations around 2104 cm^{-1} .

image pixel. The WITec Project software was used for data evaluation and processing. For further information please refer to the figure legends.




Raman images of macrophage cells incubated with 400 μM of oleic acid for different periods of time. After 30 hours macrophages store excessively lipids and develop into foam cells.

3D confocal Raman imaging of endothelial cells

In this study Raman imaging was used to three dimensionally study the heterogeneity of single endothelial cells and to define size, volume, shape and biochemical composition of the cellular organelles. The ability of confocal Raman imaging to construct 3D maps without disrupting the spatial integrity of the cell provides a unique insight into biochemical architecture and cellular processes of endothelium. 3D confocal Raman imaging can be obtained through subsequent measurement of several sample layers. And can be used e.g. for the early diagnosis of cancer by detecting subtle biochemical changes in cells and tissues associated with cancer development and progression or for morphological analysis of a tumor.

An alpha300 R confocal Raman microscope with a 60x water immersion objective for cells and an excitation wavelength of 532 nm was used for this study. The data acquisition was controlled by the WITec Project software package. All spectra were baseline corrected using a polynomial of degree 3 and the routine procedure for removal of cosmic rays was applied. The Image J processing program was applied to obtain 3D pictures of cells and tissues. For further information please refer to the figure legends.

Confocal Raman imaging-stack of an EA.hy.926 cell. Integration maps over the $\nu_{\text{C-H}}$ (2800 – 3020 cm^{-1} range).




Confocal Raman imaging-stack of an EA.hy 926 endothelial cell (A). 2D Integration maps over the CH stretching vibration (2800–3020 cm^{-1} range) at different z-positions with corresponding single spectra (B) extracted from the same point in the maps (dotted line).

Raman microscopy in malaria diagnosis


Raman microscopy can be applied as rapid and accurate diagnostic tool for malaria. Malaria is caused by the parasite Plasmodium that infects human blood cells. In the blood cells it catabolizes large amounts of haemoglobin into an insoluble bio-mineral known as haemozoin (malaria pigment). Besides serological antigen detection, this pigment is widely used for malaria diagnosis. A common microscopic method to detect haemozoin in malaria diagnosis is the darkfield microscopy. Through dark-field microscopy haemozoin is presented as bright intracellular

lar spots. The combination of darkfield and Raman microscopy enables the comprehensive characterization of these spots and provides information about the stage of the disease. In the following study Raman spectroscopy was applied to detect and image small inclusions of haemozoin in Plasmodium falciparum infected red blood cells. Therefore the confocal Raman microscope alpha300 R with a 532 nm Nd:YAG laser and a 60x/NA=1.0 water immersion objective was used. The presence or absence of haemozoin was detected by the strong 1569 cm^{-1} Raman band, which was used as marker band for haemozoin. After data acquisition the spectral data was

further processed with the WITec software. A cluster analysis was applied in order to facilitate the identification of haemozoin within the cells and to enable the classification of the different stages of disease. A cluster analysis is suitable for the automatic identification of similar spectra and the classification of multi-spectrum data into a user-defined number of clusters. Thus color-coded, user-selected clusters images can be generated. The Raman results were compared with the dark-field micrographs in order to verify the analyzing method. For more details please refer to the figure legends.



Verification of the method by comparison of dark-field and Raman microscopy (A) Visible micrograph showing the partial dark-field effect lighting up haemozoin deposits in several parasite infected cells. (B) Zoom-in of the red square in (A): The infected erythrocyte towards the bottom right corner is triply infected. (C) Chemical Raman map generated by integrating the region between 1700 and 1500 cm^{-1} . The lighter color shows regions of stronger counts. (D and E) Unsupervised Hierarchical Cluster Analysis (UHCA) map generated for the 1700 – 1500 cm^{-1} range for 2 clusters (D) and for 5 clusters (E). The haemozoin particles can be clearly identified by Raman cluster analysis.



Identification of haemozoin in different stages of disease (A) Visible partial dark-field micrograph of a thick film of malaria infected cells showing the malaria pigment haemozoin as light intracellular dots (B) Chemical map of the 1569 cm^{-1} Raman band of approximately the area bounded by the red square in (A). The yellow and red colors show the haemozoin deposits in the blood cells. (C) Cluster analysis performed using the range between 1700 cm^{-1} and 1300 cm^{-1} showing 5 clusters. The pink cluster correlates to the haemozoin deposits within late-stage of infected blood cells. The green and grey clusters are a mixture of haemoglobin and haemozoin. The light blue cluster correlates well with haemoglobin in the cell while the red-sub-micron sized dots (300 nm) appear to be inclusions of haemozoin observed in cells of different infection stages. (D) Mean extracted spectra from 5 cluster showing the region between 1600 and 1500 cm^{-1} . The colors correspond to (C) Note the differences in the red and light blue spectra. The strong shoulder on the pink spectrum centered at 1569 cm^{-1} indicates that the sub-micron-dots observed in (C) and are from inclusions of haemozoin.

Characterization of wood cells

Cell walls of wood cells consist of crystalline cellulose fibrils embedded in an amorphous hemicellulose-lignin matrix with a multilayer arrangement. The proportion of these three polymers varies between each layer. The individual layers are formed at different periods during cell differentiation. After the cell wall reaches its final structure and thickness, the mechanically crucial secondary cell wall, consisting of three different sublayers (S₁, S₂, and S₃), is formed. The cellulose microfibril arrangement and the specific chemical composition vary between sublayers. The S₂ layer is thickest

(75-85 percent of the total cell wall thickness) and most important for mechanical stability, with the highest carbohydrate content. The middle-lamella layer, mostly lignin, is attached to the primary cell wall and ensures the adhesion of a cell to its neighbouring cells. The chemical composition of the cell wall and the alignment of the cellulose microfibrils show significant interspecies and intraspecies variability. Confocal Raman microscopy was used to illustrate changes in molecular composition of secondary plant cell wall tissues of poplar (*Populus nigra* and *Populus deltoids*) wood. In this experiment, a cross section of *Populus nigra* (20 µm thick) was investigated

in water using the alpha300 R confocal Raman microscope. Figure A shows an overview of 30 x 30 µm² (150 x 150 pixels) including 22,500 spectra with an integration time of 100 ms per spectrum using a 532 nm NdYag laser for excitation. The three images in figure B-D (zoom-in at 14 x 14 µm²) were obtained by integrating over selected Raman lines. Figure E shows the color-coded image produced from the three images in figure B-D. The corresponding spectra are shown in figure F. The red spectrum represents the distribution of the G-layer, the green spectrum corresponds to the S₂-layer and the yellow spectrum is mainly found in the middle lamella.





Figure A: 30 x 30 µm² overview of *Populus nigra* in water. Figure B-D: Zoom-in, 14 x 14 µm². Integrating over different Raman lines Figure E: Color-coded image created of the three images in figure B-D. Figure F: Corresponding spectra. Red= G-layer, green= S₂-layer, yellow= middle lamella.

Imaging of banana starch

The image shows a white-light microscopy and color-coded Raman image overlay of a squashed banana pulp sample. The Raman image was acquired with the automated WITec apyrion microscope attached to the 600 mm focal length UHTS spectrometer system to provide the highest spectral performance. Red: beta-Carotene-rich areas, Green: Starch, Blue: Water. Excitation: 532 nm, Laser Power: 49.9 mW, Scan area: 400 x 300 μm^2 , 1200 x 900 pixels, Integration time: 2 ms/spectrum/pixel.



Bacteria studies

It has been reported that Raman Microscopy can be used to classify bacteria by their individual Raman spectra down to the substrain level. A high-resolution Raman imaging system also allows a single bacterium to be evaluated for metabolic products or drug detection inside the bacterium or to distinguish the intra- and inter-cellular heterogeneity. In the following study, *Legionella Bozemanii* and *Bacillus Cereus* are investigated with an



alpha300 R Confocal Raman Imaging system at the single cell level.

Legionella Bozemanii

Legionella strains can produce Poly- β -hydroxybutyric acid (PHB) in response to physiological stress, which serves as an energy storage molecule that can be detected with Raman spectroscopy. In this experiment a Raman image was acquired with a scan range of 25x25 μm^2 . The blue area in the color-coded Raman image corresponds to the integral intensity


of the C=O ester stretching band at 1726 cm^{-1} . The red area depicts the cell body imaged using the protein amide I band at 1662 cm^{-1} . The corresponding Raman spectra are shown below. The Raman image clearly reveals that the bacteria cells can contain different levels of PHB. Cells that contained little or no PHB could be discriminated from cells with a very high concentration.

(A) Raman Image of *Legionella Bozemanii* and corresponding Raman spectra (B): Vegetative cells (red) and PHB (blue).




Bacillus Cereus

In a second experiment Raman Imaging was performed on *Bacillus Cereus*. This strain is also able to produce PHB as a storage material, accumulating it as intracellular granules. (C) shows the video image indicating the region where the Raman image was acquired, which is shown in (D). The *Bacillus* strain used here is also able to form three different cellular components which can be observed in the Raman image: Vegetative cells without PHB (yellow) and vegetative cells incorporating traces of PHB (green). Additionally, *Bacillus Cereus* can form spores after entering the stationary phase of growth and accumulation of PHB. In the Raman image four spores can be seen at the right edge of the chain (magenta). The fact that all three components can be detected together highlights the diagnostic power of Confocal Raman Imaging for identifying phenotypic heterogeneity at a single cell level.




(C) Video image of *Bacillus Cereus* (D) Raman Image of *Bacillus Cereus* from the region indicated with the rectangle in (C): Vegetative cells (yellow), PHB (green), spores (magenta).

Images courtesy of Dr. Antje Hermelink, Robert-Koch-Institute, Berlin, Germany



Imaging of bacterial biofilms

In this confocal Raman study complex microbial aggregates formed of dense, highly hydrated, highly structured clusters of bacterial cells are analyzed and imaged. With confocal Raman microscopy it was possible to achieve an in-depth analysis of the targeted microbial communities without additional, more invasive techniques or prior knowledge of the sample being required. A WITec confocal Raman microscope with a 60x/NA=1.0 water immersion objective, and 532 nm excitation wavelength was used to acquire x-z-Raman depth profiles. Different bacterial micro-colonies could be distinguished and visualized. The spatial resolution was 350 nm laterally (x-direction) and 2 µm in z-direction. For further information please refer to the figure legend.



Spectral imaging based on (resonance) Raman imaging focused on the heme moiety in cytochrome-c. Set of two spectral images cutting through the outer layer of a nitrifying biofilm granule, revealing a dense multilayer of different bacterial micro-colonies: two types of *Nitrosomonas* (green and red) and one unidentified species (blue). Image dimensions: 140 x 80 spectra covering 70 x 80 µm² (x/z depth image).


Images courtesy of Dr. Ann-Kathrin Kniggendorf, Hannoversches Zentrum für Optische Technologien, Gottfried Wilhelm Leibniz Universität Hannover, Germany

Programmable Large Area Scan – Hamster brain cross section

For Raman large-area scanning on biological tissues one quite often faces a variety of obstacles. For example, fluorescence might influence the quality of the Raman spectra or the sample might not be entirely flat over the complete scan range. In order to reduce the fluorescence, a bleaching sequence before acquiring the Raman spectra sometimes eliminates the fluorescence signal. To correct for insufficient sample flatness or tilt, one approach can be to measure the sample height at selected spots. These


values can then be used to adjust the focus accordingly over the complete scan range. With the alpha300 Raman microscope, bleaching and tilt compensation can be achieved by performing a programmable large-area scan. In the following study a hamster brain cross section is moved under the laser beam and the focus position is set using the values from the predefined tilt compensation. At each image pixel a pre-bleaching sequence of 2 s is executed before the Raman spectrum acquisition which took 1 s. In (A) a video image of the complete brain cross section can be seen. (B) shows the corresponding Raman image after

evaluating the integral intensity of the CH-stretching band (a) and its center of mass (b), respectively. A second large-area scan was performed at the area indicated in (A). The results are depicted in (C) - (E) showing different images of the gray and white brain matter of that region with the corresponding spectra. This can be achieved by evaluating the Raman imaging data in more detail. Following the large-area scans a high resolution zoom-in image was acquired at the position indicated in (C), revealing differences in the chemical properties on the sub-micron scale (Fa+b).




Images courtesy of Dr. Peter Lasch, Robert-Koch-Institute, Berlin, Germany

RISE Microscopy - Correlative Raman Imaging and Scanning Electron Microscopy (Raman-SEM)



(A) Raman-SEM image overlay of a hamster brain tissue sample. In the color-coded Raman image the white brain matter is shown in Green and the gray brain matter in Red. Raman image: 100 x 100 μm^2 , 300 x 300 pixels = 90,000 spectra, 50 ms integration time per spectrum. (B) The corresponding Raman spectra reveal the different spectral characteristics of the white and gray brain matter.





Rat liver - near-field microscopy

With Scanning Near-field Optical Microscopy, it is possible to study biological subjects in liquid and in air with high resolution. Thus biological samples can be examined in their typical environments, for example in a buffer solution or special culture medium, without time consuming sample preparation. Histological microtome cuts of 130 nm thickness have been color labeled with toluidine blue. After being embedded in resin, the sample was sectioned with a microtome and floated on a cover slip. For the measurement in liquid, the cover slip was put in a petri dish with solution. The WITec SNOM objective was equipped with a sapphire window for liquid measurements. Figure A shows a 100 x 100 μm^2 overview of several liver cells (= hepatocytes, Hec) with a nucleus

(N) of 10 μm diameter. Inside the nucleus, the nucleolus is visible as a dark spot. The cells are surrounded by the intercellular clefts (I). The bright patches are hepatic sinus (S) and the little bright spots are bile canaliculi (B). Figure B is a zoom-in of a nucleus. Inside the nucleolus (Nc), euchromatin (Ec) as well as heterochromatin (Hc)


is visible. The nucleus is surrounded by the nuclear membrane (Mn). The membrane pores (Mp), which are necessary for the transfer between the cell and the nucleus, can also be seen.



(A): Hepatocyte (liver cells), overview 100 x 100 μm^2 , Hec = Hepatocyte, N = Nucleus, I = Intercellular cleft, S = Hepatic Sinus, B = Bile Canaliculi. (B) Hepatocyte nucleus, zoom-in 10 x 10 μm^2 , N = Nucleus, Nc = Nucleolus, Ec = Euchromatin, Hc = Heterochromatin, Mn = Nuclear membrane, Mp = Membrane pores.

Depth profiles of human skin

In the following study the transdermal delivery of pharmaceutical agents was investigated with Raman depth profiles of human skin tissue. Skin penetration experiments were performed by applying beta-carotene to the skin biopsies. After incubation x-z-Raman profiles of the intact full thickness skin biopsies were generated and spectral images were acquired using an alpha300 R microscope with a 50x 0.9NA objective and 785 nm excitation wavelength. Depth profiles were collected employing the x-z-scan modus starting at the bottom end of the image. The spectral information in the images highlight the different components of the tissue and the presence of the agents. Presented are Raman images of untreated skin and diffusion patterns for beta-carotene. For further information please refer to the figure legend.



Raman depth profile images of human skin. Images were collected using a 785 nm excitation at a step size of 1 µm. (A) and (B) show Raman images of normal untreated skin (B) was preprocessed using a spectral unmixing algorithm. The associated spectra plotted in (C) exhibit the typical protein characteristic peaks at 1650 – 1800 cm⁻¹. The Raman images (D) and (E) show depth profile images of skin treated with a beta-carotene formulation. The distribution of beta-carotene is plotted in red and magenta. The associated spectra plotted in (F) show the characteristic Resonance Raman bands of beta-carotene at 1524 and 1153 cm⁻¹.


Imaging of vessels with combined Raman- AFM

Vessels consist of three main parts including tunica intima, media, and adventitia. Changes within the vessel tissue structure and biochemical composition can indicate pathological alterations of the tissue, especially in endothelial cells. Those changes can be symptoms of vessel diseases e.g. diabetes, hypertension and atherosclerosis. Immunohistochemical (IHC) staining is a well-established method in pathology to analyze the cellular compartments of a sample. While for

IHC staining the sample is treated with fluorescent dyes, Raman spectroscopy distinguishes specific tissue and cell constituents by integration of the appropriate Raman band. In many cases the information obtained with these two methods is the same. In the following study, an alpha300 RA fully integrated Raman-AFM microscope was used to investigate murine vessels. The Raman-AFM cross-section analyses of the vascular wall status are useful to monitor ex vivo disease progression and effects of treatment in experimental diabetes, atherosclerosis or hypertension. While Raman

spectroscopy provides specific biochemical information about the sample composition, AFM (Atomic Force Microscopy) detects the topography, structure, and physical properties, e.g. stiffness, adhesion, of the sample's surface. Thus combined Raman- AFM analyses are useful for a comprehensive characterization of the vessel status and both pathological and physiological conditions can be investigated and distinguished. In the following figure an exemplary Raman-AFM analysis of a vessel wall cross section is shown. For further information please refer to the figure legend.

Images of the vessel wall cross section: Raman distribution image of organic vessel specimens (A). AFM AC (intermittent or tapping mode) images of the vessel topography shown in (B) and the phase image shown in (C). The phase image reveals e.g. the tissue stiffness. The 3D topography imaging (D) was used to position cross section (black and red lines). Through the cross section profiles (E) the thickness of the vessel walls can be determined. Data acquired with an alpha300 RA fully integrated microscope, 100x/NA=0.9 air objective, 532 nm excitation wavelength, 10 mW laser power at the sample. The integration time per Raman spectrum: 0.2 s. AFM topography image: AC (tapping) mode.




Detecting pathological changes in liver tissue by confocal Raman imaging

The liver is involved in most metabolic processes of the human body. Thus metabolic dysfunctions often affect the liver and cause altered tissue structures or composition. For instance a relation between diabetes mellitus type II and liver steatosis and cirrhosis has been described in literature. Furthermore liver dysfunctions may be associated with atherosclerosis. Raman spectroscopy has been widely applied to biological and biomedical samples as it has a number of useful advantages for studying

such materials: Minimal or no sample preparation, nondestructively and the possibility to gain molecular information without the use of stains. In the following study the biochemical alterations in diabetic and atherosclerotic tissue compared to normal liver tissue were investigated by confocal Raman imaging. Therefore liver tissues of atherosclerotic and diabetic mice models were used. The analyses were performed on a confocal Raman microscope alpha300 R equipped with a frequency-doubled Nd:YAG laser for 532 nm excitation and a 100x objective (NA=0.9). The Raman results were evaluated as follows: The Raman

signal of the C-H stretching vibrations ($2800 - 3020 \text{ cm}^{-1}$) originates from both lipids and proteins. To image saturated and unsaturated lipids the range between 1240 cm^{-1} and 1340 cm^{-1} was selected. The amide III region ($1215 - 1275 \text{ cm}^{-1}$) was chosen to detect proteins. Thus protein- and lipid-rich regions within the liver tissue could be determined and distinguished. It could be shown that the lipid content in the diabetic mouse tissue was highest and relatively high in atherosclerotic mice. These data strengthen the assumption that liver steatosis and metabolic diseases such as diabetes or atherosclerosis are correlated.




A microphotograph of the cross-section of liver tissues taken from control (top), diabetic (middle) and atherosclerotic (bottom) mice. The red square indicates investigated areas ($26.1 \times 21 \mu\text{m}^2$, $18.8 \times 15.5 \mu\text{m}^2$ and $19.2 \times 16.7 \mu\text{m}^2$ respectively). Raman maps of the lipid and protein distributions obtained by integration of the marker bands in the regions of: $3020-2800 \text{ cm}^{-1}$ (lipids and proteins), $1240-1340 \text{ cm}^{-1}$ (lipids) and $1215-1275 \text{ cm}^{-1}$ (amide III, proteins). Scale bars as presented in each image.


Histopathological human breast cancer tissue

For histopathological diseases or cancer detection, Raman imaging is considered to be a potential alternative for the currently used techniques, which require staining of the biopsy tissue before the manual microscopic investigation. As Raman imaging is a label free characterization method for the sample's compounds, it is well-suited to generally accelerate the detection process or to establish automated and more reliable detection procedures. However, extended basic research must still be made in order to achieve clinical relevance for Raman imaging as a standard method for diagnosis. In the following study untreated histopathological breast cancer tissue was investigated with the WITec alpha500 Confocal Raman Microscope using its large area scanning capabilities. (A) shows an overview video image of the sample. Based on this image, the scan ranges for the Raman imaging scans can be selected as indicated by the rectangles. The first large-area scan (red rectangle) was acquired in order to obtain an overview Raman image of the tissue sample followed by a second large area scan (green rectangle) at the area with potentially carcinogenic tissue. (B) - (D) show the resulting Raman Images and corresponding spectra. The imaging parameters for the first large area scan are 10.2 x 9.7 mm², 150 x 150 pixels (= 22,500 spectra). The scan range of the second scan was 4.85 x 4.00 mm² and 150 x 150 pixels. Focusing on (D) reveals differences in the chemical composition of the tissue. Especially the tissue in the areal indicated with the orange area shows clear differences in the spectral characteristics (no red spectrum detectable). Using staining techniques (results not shown here) it could be shown that the carcinogenic


tissue is located exactly at this area of the image. This results suggest, that a differentiation of healthy and carcinogenic tissue with Raman Imaging can be possible.




(A) Video Image of the unstained breast cancer tissue area indicates the location of the carcinogenic tissue.



(B) Large-area scan I (red rectangle in A)



(C) Basis spectra used for image generation for the images in (B) + (D)









(D) Large-area scan II (green rectangle in A). The orange area indicates the location of the carcinogenic tissue.

Comparison of noncancerous and cancerous tissue

Through confocal Raman imaging the general tissue composition as well as the appearance of tumor markers can be identified. Thus noncancerous and cancerous tissue can be distinguished and characterized without any external agent and stain or sample prepara-

tion prior to the experiments. In this study noncancerous and cancerous human breast tissue of the same patient was examined with an alpha300 R confocal Raman microscope, a 50x objective, and a 532nm excitation laser with 10 mW laser intensity at the sample. The WITec Project plus software was used for data acquisition and processing. The appearance

of different tissue components was determined through their unique Raman spectra and the noncancerous and cancerous tissues could be distinguished by the altered appearance of carotenoids, lipids, and proteins. For further information please refer to the figure legends.









carotenoids

lipids

proteins

fluorescence



carotenoids

lipids

proteins

fluorescence

Breast tissue from the margin of the tumor mass:

Raman image ($350 \times 350 \text{ mm}^2$), integration time 0.03 sec, 1 accumulation; average spectra used for the basis analysis method, colors of the spectra correspond to the colors of the different areas in the Raman image; images for the filters for spectral regions: carotenoids ($1490 - 1580 \text{ cm}^{-1}$), lipids ($2850 - 2950 \text{ cm}^{-1}$), proteins ($2900 - 3010 \text{ cm}^{-1}$), and fluorescence ($2200 - 2300 \text{ cm}^{-1}$). Please note that carotenoids are clearly detectable in the tissue.

Breast tissue from the tumor mass (carcinoma ductale G3 infiltrans mammae): Raman image ($150 \times 150 \text{ mm}^2$), integration time 0.03 sec, 1 accumulation; average spectra used for the basis analysis method, colors of the spectra correspond to the colors of the different areas in the Raman image; images for the filters for spectral regions: carotenoids ($1490 - 1580 \text{ cm}^{-1}$), lipids ($2850 - 2950 \text{ cm}^{-1}$), proteins ($2900 - 3010 \text{ cm}^{-1}$), and fluorescence ($2200 - 2300 \text{ cm}^{-1}$). Please note that the amount of detectable carotenoids is strongly reduced.

Images courtesy of Halina Abramczyk and Beata Brozek-Pluska, Lodz University of Technology, Institute of Applied Radiation Chemistry, Laboratory of Laser Molecular Spectroscopy, Wroblewskiego 15, 93-590 Lodz, Poland.