

## UHP-T-EP

### Ultra-High Power LED Light Source

Ver. 03

#### Introduction

The Ultra High Power LED provides more than X10 increase of power in comparison with similar High power LED devices. The Ultra High Power LED (> 50 **Watt**) is an effective replacement of lasers and lamps in many power demanding applications, such as fluorescence microscopy, Optogenetics in-vitro, patterned illumination, confocal microscopy, chemical reaction activation and numerous others. Fluorescence microscopy applications will benefit from highly homogenous and flat field illumination. The shielded UHP-T-EP LED head contains the high current driver, while the LED controller box contains control functions. This arrangement eliminates much of RFI / EMI interference common in high current light sources. The controller features Optically Isolated TTL and Analog Inputs. These features make this product especially suited for electrophysiology rig applications. UHP-T-EP can be used with a full range of other optical OptiBlocks such as Fiber Coupler, Liquid Light Guide coupler, Beam Combiner, Filter Wheel, Beam Switcher and others. See Optional Accessories section below.



#### Key Features

- Single ultra-high brightness LED chip provides highly homogeneous illumination over whole field of view.
- Shielded LED head with high current driver, low RFI/EMI.
- Optically isolated TTL input for external triggering (no shutter needed) or strobe operation
- Optically isolated Analog input (0-5V) for LED power control by external device like D/A interface
- Computer control via USB by Windows software, LabView VI or uManager (optional).
- Long life (no lamp or laser tube replacement required)
- Compatible with Prizmatix modular Microscope-LED Light Source products family – see below for details.
- Special low optical noise model available for detection of small signals (option).
- Stable precisely adjustable power

---

#### Elliot Scientific Limited

Unit 11 Sandridge Park, Porters Wood, St Albans, AL3 6PH, United Kingdom  
Tel: +44 (0)1582 766300 Fax: +44 (0)1582 766340  
[www.elliotscientific.com](http://www.elliotscientific.com) [sales@elliotscientific.com](mailto:sales@elliotscientific.com)


# Prizmatix

## Applications

- Fluorescence microscopy in electrophysiology rig.
- *In-Vitro* Optogenetics
- Multi-wavelength systems
- Bio analysis
- OEM

## Optical Specifications – Color LEDs

| P/N          | Peak<br>[nm] | Centroid<br>[nm] | FWHM<br>[nm] | Power Collimated<br>[mW] |
|--------------|--------------|------------------|--------------|--------------------------|
| UHP-T-405-EP | 404          | 406              | 16           | 2100                     |
| UHP-T-450-EP | 451          | 453              | 25           | 2300                     |
| UHP-T-520-EP | 520          | 527              | 37           | 1000                     |
| UHP-T-625-EP | 624          | 620              | 24           | 850                      |
| UHP-T-630-EP | 636          | 632              | 26           | 1000                     |
| UHP-T-650-EP | 656          | 652              | 26           | 1000                     |



### Elliot Scientific Limited


Unit 11 Sandridge Park, Porters Wood, St Albans, AL3 6PH, United Kingdom

Tel: +44 (0)1582 766300 Fax: +44 (0)1582 766340

[www.elliotscientific.com](http://www.elliotscientific.com) [sales@elliotscientific.com](mailto:sales@elliotscientific.com)



# Prizmatix



## Benchtop LED Current Controller Specifications

- Constant current or chopping modes
- Precise LED current setting by 10 turn dial
- TTL external trigger input
- Analog input for external LED power control (0-5Vdc)
- Optically isolated TTL and Analog inputs
- Compact and robust enclosure



|                                     |         |                           |
|-------------------------------------|---------|---------------------------|
| Connectors for TTL and Analog input |         | Optically Isolated BNC    |
| Digital modulation frequency        | Hz      | DC - 10000                |
| Rise / Fall time                    | $\mu$ s | <10 / <3                  |
| Analog input voltage range          | V       | 0 - 5                     |
| Analog modulation frequency         | Hz      | DC - ~100                 |
| Current controller supply voltage   | V       | 12                        |
| Power adaptor input                 |         | 85-264 VAC, 47-63Hz, 1.5A |

## General Specifications

|                                              |    |          |
|----------------------------------------------|----|----------|
| Operation temperature range                  | °C | 10 - 35  |
| Storage temperature range                    | °C | -10 - 55 |
| Operating relative humidity (Non condensing) | %  | <90      |

### Elliot Scientific Limited


Unit 11 Sandridge Park, Porters Wood, St Albans, AL3 6PH, United Kingdom  
Tel: +44 (0)1582 766300 Fax: +44 (0)1582 766340  
www.elliotscientific.com sales@elliotscientific.com



# Prizmatix

|                                      |     |                                                                                    |
|--------------------------------------|-----|------------------------------------------------------------------------------------|
| Head dimensions                      |     | See drawing below                                                                  |
| Head weight                          | g   | 350                                                                                |
| Controller dimensions (L x W x H)    | mm  | 197 x 174 x 80                                                                     |
| Controller weight                    | g   | 400                                                                                |
| Power adaptor dimensions (L x W x H) | mm  | 175 x 72 x 35                                                                      |
| Power adaptor weight                 | g   | 650                                                                                |
| Power Adaptor Safety                 |     |  |
| LED Head fan noise                   | dBA | 38                                                                                 |

## Mechanical Drawings



\* Specifications subject to changes without notice

### Elliot Scientific Limited

Unit 11 Sandridge Park, Porters Wood, St Albans, AL3 6PH, United Kingdom

Tel: +44 (0)1582 766300 Fax: +44 (0)1582 766340

www.elliotscientific.com sales@elliotscientific.com




# Prizmatix

## Optional Accessories

For full details on optional accessories please see:

<http://www.prizmatix.com/optogenetics/Optogenetics-LED-Light-Sources-and-Fiber-Optics.htm>



### Beam Combiner [2]:

Multiple LED beams can be combined into one output beam. For example UV LED can be combined with White LED to create Mercury lamp like configuration. For more details please see:

<http://www.youtube.com/watch?v=iv7dlwLHaUE>

### Filter Wheel [7]:

The UHP-Mic-LED can be equipped with a 6 positions filter wheel at the beam output. This accessory is especially useful for UHP-Mic-LED-White light source.

Please see <http://www.prizmatix.com/Optics/filter-wheel.htm> for more details.

### Fiber Coupler Adaptor [3]:

The UHP-Mic-LED can be easily changed from direct microscope coupling to fiber coupled LED configuration by means of Fiber Coupler Adaptor (SMA, CF or ST connector). Please see video clip <http://www.youtube.com/watch?v=iv7dlwLHaUE> for more details.

### Liquid Light Guide Adaptor [9]:

The Microscope-LED can be easily changed from direct microscope coupling to Liquid Light Guide coupled LED configuration by means of LLGA Adaptor. Please see video clip <http://www.youtube.com/watch?v=iv7dlwLHaUE> for more details.

### Fiber Optics Collimator [13]:

The output from optical fiber is divergent according to fiber NA. In order to reduce the divergence angle a collimator module can be used. Prizmatix collimator was especially designed to fit thick core high NA Polymer Optical Fibers. See more details at: <http://www.prizmatix.com/Optics/collimator.htm>

### Fiber Bundles [4]:

To combine outputs of multiple LEDs a Y-shaped fiber bundle with two or more input branches can be used. Prizmatix can help to configure and build custom fiber bundles for specific applications. See more details at: <http://www.prizmatix.com/dev/Custom-Fiber-Optic-Assemblies.htm>

## Elliot Scientific Limited

Unit 11 Sandridge Park, Porters Wood, St Albans, AL3 6PH, United Kingdom  
Tel: +44 (0)1582 766300 Fax: +44 (0)1582 766340  
[www.elliotscientific.com](http://www.elliotscientific.com) [sales@elliotscientific.com](mailto:sales@elliotscientific.com)

