





## GRANULE SECRETION - EXOCYTOSIS




Insulin secretion is critical to regulate blood sugar levels and any dysregulation can have significant implications for disease. Upon detection of high blood sugar insulin is secreted, stimulating cells to absorb sugar from the blood, effectively lowering the blood sugar level. Insulin is stored within granules that fuse with the plasma membrane upon stimulation. This important biological process is known as exocytosis.

In the image to the left, the granule content has been fluorescently labeled. Upon stimulation with glucose the cargo is released, thereby indirectly visualizing the release of insulin. This first appears as a faint red foci that intensifies as the vesicle approaches the cell membrane. As the granule fuses with the membrane the content leaks out of the cell, resembling a firework-like event.



The image is a still frame taken from the timelapse available on our webpage: [chipnano.com](http://chipnano.com). The granules fill the entire cell volume making it necessary to use total internal reflection fluorescence microscopy (TIRFM) to ensure that only the fluorescent signal at the cell surface is visible. Additionally, due to the heterogeneity of cells and the downstream signaling effects of insulin, the release-per-cell itself can vary greatly. ACP ZERO achieves the largest TIRF field of view. The image above achieves a uniformly illuminated field of view across 260 x 260 μm, with a time resolution of 30 frames per second. Below is a magnified timelapse of a granule secretion event.



**ACP ZERO STANDARD**

**TECHNICAL SPECIFICATIONS**

**WORKING PRINCIPLE**

**ACP ZERO** is an accessible and robust TIRF imaging platform. Our multimode waveguide technology decouples the illumination and imaging pathways. The imaging path follows conventional standards, while the illumination path occurs via built-in waveguides on the surface of a photonic chip. There are several parallel waveguides on top of the chip, allowing for imaging within an active waveguide while preserving the neighboring fluorophores. The sample is placed or grown directly on top of the chip. After automatic laser coupling into a selected waveguide, an evanescent field penetrates the sample and excites the fluorophores which are then detected by an upright microscope.

**SYSTEM SPECIFICATIONS**

Standard model, customizable options

*Camera unit*

Sensor type: Scientific CMOS  
Pixel count: 3200x3200

*Laser unit*

Standard wavelengths: 491 nm, 561 nm, 640 nm  
Customizable laser lines: TBD

*White light source (EPI)*


Technology: 200W metal halide

*Objective lenses*

Lens turret: 2 objectives, simultaneously mounted  
Lens selection:  
- 60x NA 1.2W  
- 25x NA 0.8  
- 20x NA 0.45  
- 10x NA 0.4

Lateral resolution: 50 nm @ 60x (SMLM)  
75 nm @ 25x (SMLM)  
140 nm (SMLM) @20x  
Diffraction limited and DV enhanced

Software: Acquisition software  
Post-processing suite with GPU-accelerated visualization and reconstruction for SMLM data, and deconvolution software



**CONTACT INFORMATION**

Chip NanoImaging AS  
Visiting address:  
Fiolvegen 15  
9016 Tromsø  
Norway  
cniptest@chipnano.com  
chipnano.com/contact/

chipnano is a trademark of Chip NanoImaging AS.  
© Chip NanoImaging AS.  
November 2024