


The MadMotor™ UHV can be used for coarse positioning of UHV nanopositioners such as MCL's Nano-HS3

4xC-Bored for M2.5 SHCAP on a 54mm x 54 mm square



MadMotor™ UHV

Dimensions are in inches[mm]
 Drawing No.
 MMUHV2.14.A

DRAWING IS TO BE USED AS A GUIDE ONLY. DIMENSIONS ARE SUBJECT TO CHANGE WITHOUT NOTICE. THE INFORMATION CONTAINED IN THIS DRAWING IS THE SOLE PROPERTY OF MAD CITY LABS. ANY REPRODUCTION IN PART OR AS A WHOLE IS PROHIBITED.