

## OPTOSci Price List

### Laser Current & TEC Driver Modules

*Valid from January 2016*

Models	Price	Specifications
<b>LDR250GKE</b>	<b>£ 759</b>	Eurocard, grounded cathode device current and TEC controller module for demountable 14-pin Butterfly package with USB interface, 250 mA diode drive and 2.5 A TEC drive.
<b>LDR250GAE</b>	<b>£ 771</b>	Eurocard, grounded anode device current and TEC controller module for demountable 14-pin Butterfly package with USB interface, 250 mA diode drive and 2.5 A TEC drive.
<b>LDR1000E</b>	<b>£ 771</b>	Eurocard device current and TEC controller module for demountable 14-pin butterfly package with USB interface, 1 A diode drive and 2.5 A TEC drive.
<b>LDR1500E</b>	<b>£ 820</b>	Eurocard device current and TEC controller module for demountable 14-pin butterfly package with USB interface, 1.5 A diode drive and 4 A TEC drive.
<i>The above units are supplied with full V-Drive control software, USB interface, PC communications cable, appropriate mains power supply, operators manual, heatsink and ZIF mount for 14-pin butterfly packaged fibre coupled laser, SOA or similar device. For other options and device mounts please contact us.</i>		
<b>Option</b>		
<b>-IF</b>	<b>£ 68</b>	External InGaAs photodiode in FC receptacle integrated on board.

PRICES ARE EX WORKS HARPENDEN. SHIPPING CHARGE AND VAT EXTRA  
PRICES SUBJECT TO CHANGE WITHOUT NOTICE

#### **Elliot Scientific Ltd**

3 Allied Business Centre, Coldharbour Lane, Harpenden, Hertfordshire, AL5 4UT, UK  
Tel: +44 (0)1582 766 300 Fax: +44 (0)1582 766 340 Eml: sales@elliotscientific.com