





REVISONS			
REV.	DESCRIPTION	DATE	APPROVED



3° ROTATION AXES  
±1mm X&Y TRANSLATION,  
±5mm Z TRANSLATION



GENERAL VIEW SHOWING ROTATION  
AND TRANSLATION AXES

PROPRIETARY AND CONFIDENTIAL  
COPYRIGHT ELLIOT SCIENTIFIC Ltd.

THE INFORMATION CONTAINED IN THIS  
DRAWING IS THE SOLE PROPERTY OF  
ELLIOT SCIENTIFIC LTD. REPRODUCTION  
IN PART OR AS A WHOLE WITHOUT WRITTEN  
PERMISSION IS PROHIBITED.

DIMENSIONS ARE IN mm  
GENERAL TOLERANCES: ± 0.1  
ANGULAR TOLERANCES: ±  
SURFACE FINISH:  
ALL BURRS, SHARP EDGES  
AND CORNERS TO BE  
REMOVED

AUTHOR	NAME	DATE
CHECKED	GW	01/06/2010
MATERIAL	---	
FINISH	---	
DO NOT SCALE DRAWING		

**Elliot Scientific**

TITLE  
**5 AXIS MICROPOSITIONER**

SIZE **A4** DWG. NO. **MDE277M**

SCALE:1:1 THIRD ANGLE PROJECTION SHEET 1 OF 1