




REV.	REVISIONS	DESCRIPTION	DATE	APPROVED



2 HOLES FOR MOUNTING TO FLEXURE STAGE



PLATFORM SUPPLIED WITH 2 MOUNTING SCREWS, CLAMP SET AND HEX KEY

**PROPRIETARY AND CONFIDENTIAL**  
**COPYRIGHT ELLIOT SCIENTIFIC Ltd.**  
 THE INFORMATION CONTAINED IN THIS DRAWING IS THE SOLE PROPERTY OF ELLIOT SCIENTIFIC LTD. REPRODUCTION IN PART OR AS A WHOLE WITHOUT WRITTEN PERMISSION IS PROHIBITED.

DIMENSIONS ARE IN mm  
 GENERAL TOLERANCES:  $\pm 0.1$   
 ANGULAR TOLERANCES:  $\pm$   
 SURFACE FINISH:  
 ALL BURRS, SHARP EDGES AND CORNERS TO BE REMOVED

NAME	DATE
AUTHOR GW	12/07/2010
CHECKED --	--

MATERIAL  
ALUMINIUM ALLOY

FINISH  
ANODISED BLACK

DO NOT SCALE DRAWING

**Elliot Scientific**

TITLE  
**FIXED PLATFORM, LARGE**

SIZE **A4** DWG. NO. **MDE147E**

SCALE:1:1 THIRD ANGLE PROJECTION SHEET 1 OF 1