



EACH AXIS HAS  $\pm 2.5\text{mm}$  OF TRAVEL




GENERAL VIEW  
SCALE: 1:1

	REVISIONS			
REV.	DESCRIPTION		DATE	APPROVED

PROPRIETARY AND CONFIDENTIAL  
COPYRIGHT ELLIOT SCIENTIFIC Ltd.  
THE INFORMATION CONTAINED IN THIS  
DRAWING IS THE SOLE PROPERTY OF  
ELLIOT SCIENTIFIC LTD. REPRODUCTION  
IN PART OR AS A WHOLE WITHOUT WRITTEN  
PERMISSION IS PROHIBITED.

DIMENSIONS ARE IN mm  
GENERAL TOLERANCES:  $\pm 0.1$   
ANGULAR TOLERANCES:  $\pm$   
SURFACE FINISH:  
ALL BURRS, SHARP EDGES  
AND CORNERS TO BE  
REMOVED

		NAME	DATE	<div></div>	
AUTHOR		GW	11/05/2010		
CHECKED		--	--		
MATERIAL STAINLESS STEEL, ALUM. ALLOY				TITLE  2 AXIS DOVETAIL SLIDE	
FINISH --				SIZE <b>A4</b>	DWG. NO. <b>MDE262</b>
DO NOT SCALE DRAWING				SCALE:1:1	THIRD ANGLE PROJECTION
				SHEET 1 OF 1	