






REV. NO.		REVISIONS	DATE	APPROVED
1		DESCRIPTION		
2				

4 HOLES M2 x 2 DEEP



GENERAL VIEW
SCALE: 1:1



EACH AXIS HAS ±2.5mm OF TRAVEL

0.5mm PITCH MICROMETERS
50 DIVISIONS PER TURN

PROPRIETARY AND CONFIDENTIAL COPYRIGHT ELLIOT SCIENTIFIC Ltd. THE INFORMATION CONTAINED IN THIS DRAWING IS THE SOLE PROPERTY OF ELLIOT SCIENTIFIC LTD. REPRODUCTION IN PART OR AS A WHOLE WITHOUT WRITTEN PERMISSION IS PROHIBITED.		DIMENSIONS ARE IN mm GENERAL TOLERANCES: ± 0.1 ANGULAR TOLERANCES: ± SURFACE FINISH: ALL BURRS, SHARP EDGES AND CORNERS TO BE REMOVED		Elliot Scientific	
AUTHOR GW 14/02/2007		NAME GW 14/02/2007		DATE 14/02/2007	
CHECKED --		MATERIAL STAINLESS STEEL, ALUM. ALLOY		TITLE 2 AXIS DOVETAIL SLIDE	
FINISH --		SIZE A4		DWG. NO. MDE262M	
DO NOT SCALE DRAWING		SCALE:1:1		THIRD ANGLE PROJECTION	
				SHEET 1 OF 1	