




REV. /		REVISIONS	DATE	APPROVED
REV.		DESCRIPTION		



GENERAL VIEW
SCALE 2:1



<div>PROPRIETARY AND CONFIDENTIAL</div> <div>COPYRIGHT ELLIOT SCIENTIFIC Ltd.</div> <div>THE INFORMATION CONTAINED IN THIS DRAWING IS THE SOLE PROPERTY OF ELLIOT SCIENTIFIC LTD. REPRODUCTION IN PART OR AS A WHOLE WITHOUT WRITTEN PERMISSION IS PROHIBITED.</div>	DIMENSIONS ARE IN mm		MATERIAL		TITLE	
	GENERAL TOLERANCES: ± 0.1		ALUMINIUM ALLOY		ROTATION STAGE ADAPTER	
	ANGULAR TOLERANCES: ±		FINISH			
	SURFACE FINISH:		ANODISED CLEAR		SIZE	DWG. NO.
	ALL BURRS, SHARP EDGES AND CORNERS TO BE REMOVED		DO NOT SCALE DRAWING		A4	MDE859
				SCALE:2:1	THIRD ANGLE PROJECTION	SHEET 1 OF 1